



ITEM NO.	PART NUMBER	DESCRIPTION	Qty.
1	HT-R12-ELE202	PLATE, ELEC. BOX TOP	1
2	HT-R12-ELE001	PLATE, ELECTRICAL BOX SIDE	1
3	HT-R12-ELE009	PLATE, ELEC. BOX FRONT PNL	1
4	HT-R12-ELE010	GROMMET, ELEC. BOX	1
5	*SEE NOTES*	TRANSFORMER	1
6	HT-R12-RLY	RELAY	1
7	HT-R12-ELE100	PLATE, COMPONENT MTG.	1
8	HT-R12-HLW	HTR/ HIGH LIMIT WASHER	1
9	*SEE NOTES*	HEATING ELEMENT	1
10	** MPS-30-12**	12V DC POWER SUPPLY	1
11	**HT-R12-ELE102**	BRKT, POWER SUPPLY	1

NOTES:
 MAINS VOLTAGE DEPENDANT 115V-120V, 220-240V
 HEATING ELEMENT: HT-R12-120HE (120V)
 HEATING ELEMENT: HT-R12-220HE (220V)
 TRANSFORMER 115-120V: HT-R12-120XFM
 TRANSFORMER 220-240V: HT-R12-220XFM

** DELUXE MODELS ONLY**

UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES TOLERANCES: ANGLES: ±1 FRACTIONS: ±1/16 DECIMALS: .XX ±.03 .XXX ±.005		
DRAWN	NAME TPH	DATE 10/23/2016
ENGINEERING APPR.		
Q.A.		
PROPRIETARY AND CONFIDENTIAL THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF RICHMAR-NAIMCO. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF RICHMAR-NAIMCO IS PROHIBITED.		
MATERIAL		
FINISH		
DO NOT SCALE DRAWING		



TITLE:
 ELECTRICAL ENCLOSURE ASSY

SIZE A	DWG. NO.	REV B
------------------	----------	----------

NOTED SHEET 5 OF 7