



ITEM NO.	PART NUMBER	DESCRIPTION	QTY.
1	HT-R12-TUBTHM	THERMOSET TUB	1
2	HT-R12-HBF	3/8 HOSE BARB FITTING	1
3	HT-R12-HLN	HIGH LIMIT NUT	1
4	HT-R12-HLW	HTR/ HIGH LIMIT WASHER	1
5	HT-R12-HBW	3/8" HOSE BARB FITTING WASHER	1
6	HW-520BRNUT	NUT, 1/2" HOSE BARB	1
7	HW-HCLP	HOSE CLAMP W/1/4" NUT	2
8	HT-R12-TUB-001	BRKT, TUB TO BASE PLATE MTG	1
9	HT-R12-SEAL	GASKET, UPPER HYDROTHERM	1
10	HW-CHOSE	DRAINLINE	6 FT.
11	HT-R12-HLF	HIGH LIMIT FITTING	1
12	HT-R12-HLS	HIGH LIMIT SWITCH W/HARNNESS	1
13	HW-DHV	DRAIN HOSE VALVE	1

UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES TOLERANCES: ANGLES: ±1 FRACTIONS: ±1/16 DECIMALS: .XX ±.03 .XXX ±.005		
DRAWN	NAME	DATE
	RA	9/12/2016
ENGINEERING APPR.		
Q.A.		
PROPRIETARY AND CONFIDENTIAL THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF RICHMAR-NAIMCO. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF RICHMAR-NAIMCO IS PROHIBITED.		
MATERIAL		
FINISH		
DO NOT SCALE DRAWING		

TUB ASSEMBLY

SIZE	DWG. NO.	REV
A		B
NOTED		SHEET 4 OF 7

5

4

3

2

1